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PAG
3.8

Sample Prep Solution MOSH/MOAH

MOSH/MOAH Sample Prep Solution performs fully automated sample preparation and sample introduction for efficient determination of mineral oil residues in extracts of food, feed, body care products and packaging. The system is based on an online-coupled HPLC-GC/FID system using for automated sample preparation and introduction.

In the initial LC step, mineral oil residue is separated into two fractions: Mineral oil saturated hydrocarbons (MOSH) and mineral oil aromatic hydrocarbons (MOAH). These fractions are subsequently transferred to two separate GC columns for individual analysis in a combined dual channel GC system.

1 Online HPLC-GC/FID system

The MOSH/MOAH Sample Prep Solution is based on an HPLC system coupled directly to a dual channel GC/FID system. The configuration shown here includes optional hardware for the Epoxidation and AlOx cleanup steps.

2 MultiPurpose Sampler (MPS)

The MPS is used for sample preparation and -introduction and can also perform fraction collection.

The solution meets the requirements of the DIN EN 16995:2017-08 Standard. The dual channel GC separation with FID detection enables a complete MOSH/MOAH analysis in only 30 minutes.

LC- and GC chromatograms are displayed in real time to facilitate method development or method optimization, for example, adjusting switching times for MOSH and MOAH fraction collection.

No additional control box is required to set up and control method parameters such as carrier gas pressure or flow.

The combined software controls all method parameters for sample preparation, HPLC and GC. The entire workflow is conveniently and efficiently processed under one integrated user interface and sequence table.

The MOSH/MOAH Sample Prep Solution can be extended to perform additional widely used sample preparation steps depending on your requirements:

- Epoxidation to remove interfering naturally occurring olefins
- ALOX Clean-up to retain and remove long chain n-alkanes of plant origin

3 Fractionation by HPLC

The MOSH and MOAH fractions are separated on an HPLC column, which retains triglycerides. During the ensuing GC run, the HPLC column is back-flushed to remove the triglycerides. When the GC is ready, the HPLC column is clean and ready for the next sample.

MOSH/MOAH Sample Prep Solution delivers

Fast answers and reliable results based on efficient automation and intuitively operated software.


- Determination of MOSH and MOAH in 30 minutes
- Efficient processing and reporting of MOSH/MOAH data, including data generated with other analysis systems.
- Specially developed algorithms for hump- and peak detection ensure fast and efficient data processing and correct results
- Manual method adjustments and reintegration are easily performed at any time.

Simple and efficient operation

- Method modification by mouse-click
- Unique transverse mounting of GC columns enables easy access and simplified maintenance
- Integrated user interface for the complete system

Flexible adaptation to individual requirements

- Modular addition of Sample Preparation steps, such as, for example, automated epoxidation or ALOX clean-up
- Easy customization of report formats

GC oven

Perfectly organized: In the GC oven, the retention gaps (A) and GC columns (B) are transverse mounted and positioned for easy access.

PAG
6.8


Analysis details

The MOSH and MOAH fractions are separated by HPLC and individually transferred to their designated GC channel. The HPLC separation is performed using a normal phase silica gel column and n-hexane/dichloromethane mobile phase.

The individual 450 µL fractions are transferred to their respective GC channels while triglycerides are retained on the HPLC column.

During the GC run, the HPLC column is back-flushed and cleaned. This ensures that correct results can be obtained for the following sample while safeguarding reliable system operation.

In the Early Vapor Exit, excess HPLC eluent is removed before the MOSH and MOAH fractions are transferred to their respective GC separation columns.

PAE
7.8


MOSH/MOAH analysis results in 4 easy steps

The MOSH/MOAH Data Processing Software uses dedicated algorithms developed specifically for hump- and peak detection

to ensure efficient data processing as well as accurate results and reporting. Manual adjustments and batch reintegration are easily performed at any time.


Data import 1

Batch processing of a large number of chromatograms using standardized conditions ensures high throughput

2 Automated sample type recognition


Chromatogram file information helps determine whether files represent samples or standards and is used for sample description

3

Visual inspection of integration intervals

Predefined integration intervals for standard alkanes are clearly displayed and can be adjusted as required, enabling fast visual inspection and reliable identification of standard peaks

4 Integration interval selection

Integration windows for MOSH and MOAH and for standards are selected by mouse-click

Name	MOSH	MOAH	Start (Minutes)	Center (Minutes)	Stop (Minutes)
△ C11	<input type="checkbox"/>	<input type="checkbox"/>	11.209	11.422	11.634
△ C13	<input type="checkbox"/>	<input type="checkbox"/>	12.396	12.614	12.832
△ C15	<input type="checkbox"/>	<input type="checkbox"/>	14.407	14.515	14.625
△ C17	<input type="checkbox"/>	<input type="checkbox"/>	16.529	16.579	16.628
△ C19	<input type="checkbox"/>	<input type="checkbox"/>	18.290	18.332	18.373
△ C21	<input type="checkbox"/>	<input type="checkbox"/>	19.687	19.729	19.772
△ C23	<input type="checkbox"/>	<input type="checkbox"/>	19.759	19.847	19.935
△ C25	<input type="checkbox"/>	<input type="checkbox"/>	21.909	22.059	22.218
△ C27	<input type="checkbox"/>	<input type="checkbox"/>	23.303	23.392	23.481
△ C29	<input type="checkbox"/>	<input type="checkbox"/>	25.491	25.590	25.629